




Heat and Moisture Exchangers Moistrap



- Installed Geared Swivel Connector.
- L & M, 2 sizes, Straight & Bent type.
Every HME has sampling port.



Geared swivel connector

Features

Moistrap M-S3.0/M-S3.5

- HME directly connected with Tracheal tube.
- Suction port is lined up straightly with connector of Tracheal tube

Moistrap M-S15

- HME applicable for the patients under inhalation Anesthesia system or Respirator.
- Gas monitoring is available through sampling port.



Moistrap M-MS / M-MB / M-LS / M-LB

- HME installed Geared Swivel Connector.
- Corresponded to the ventilation volume, 2 sizes and straight & bent type are available. Every HME has sampling port.



Moistrap M-MS
(straight type)



Moistrap M-MB
(bent type)



Moistrap M-LS
(straight type)



Moistrap M-LB
(bent type)

Specification

Spec./Model	M-S3.0	M-S3.5	M-S15	M-MS	M-MB	M-LS	M-LB
Tidal volume (mL)	10~50	10~50	15~75	50~600	50~600	150~1500	150~1500
Dead space (mL)	2.7	2.8	3.1	11	11	25	24
Moisture loss (mg±0.1air)	7.6 (Vt 50mL)	7.6 (Vt 50mL)	9.4 (Vt 70mL)	10.9 (Vt 500mL)	10.9 (Vt 500mL)	7.9 (Vt 500mL)	7.9 (Vt 500mL)
Pressure drop (hPa)	1.9 (5L/min)	1.4 (5L/min)	I 0.7 (10L/min)	I 0.8 (30L/min)	I 0.9 (30L/min)	I 0.5 (30L/min)	I 0.5 (30L/min)
			E 0.7 (10L/min)	E 0.8 (30L/min)	E 0.9 (30L/min)	E 0.4 (30L/min)	E 0.5 (30L/min)
Weight (g)	4	4	6	13	14	23	24
Connector size (mm)	Patient side OD3.0	Patient side OD3.5	Patient side ID15	Patient side ID15 OD22	Patient side ID15 OD22	Patient side ID15 OD22	Patient side ID15 OD22
	Machine side OD15	Machine side OD15	Machine side OD15	Machine side OD15	Machine side OD15	Machine side OD15	Machine side OD15
Sterilization	EOG	EOG	EOG	EOG	EOG	EOG	EOG
Package (box)	20pcs	20pcs	20個	20pcs	20pcs	20pcs	20pcs