



Heat and Moisture Exchangers **Sofitvent**



- Applicable from Infant to Adult.
- Will support more active motion of the patients.

Features



Sofitvent SV-L with Sofit clear NC

- HME connected with Tracheostomy tube for the patients from Infant to Adult.
- It is easy to handle when you hold the dent of the body.



Sofitvent SV-LO2 with Sofit clear NC and oxygen tube

- Sofitvent SV-LO2 is available to supply oxygen to the patients who need oxygen therapy.
- It is available to suck the sputum by suction catheter through the slit.



Sofitvent SV-S with Sofit P

- By low height design, Sofitvent SV-S is available for neonates and infants.
- Light weight reduces the load on the patient.

Specification

Spec./Model	SV-L	SV-LO2	SV-S
Tidal volume (mL)	50~	50~	35~250
Dead space (mL)	10	10	6
Moisture loss (mgH ₂ O/L air)	12.9(Vt 500mL)	12.9(Vt 500mL)	13.6(Vt 250mL)
Pressure drop (hPa)	0.27(30L/min)	0.27(30L/min)	0.16(10L/min)
Weight (g)	5	6	3
Connector size (mm)	Patient side ID 15	Patient side ID 15	Patient side ID 15
Sterilization	EOG	EOG	EOG
Package (box)	30pcs	30pcs	30pcs