



Sterile Water for Inhalation Calm Pure

Calm Pure (sterile water bag) is closed humidifier system for oxygen inhalation which will be contributed to efficient nursing work and infection control in the hospital.

- Calm Pure can be used as humidifier (oxygen therapy) or nebulizer(aerosol therapy) by selecting 2 different sterilized adapters with same bag.
- Water bag pre-filled with endotoxin-free sterile water is sterilized by autoclave.
- Capacity of bag, 250ml, 500ml, 800ml, 3 sizes are available. Size will be selected in accordance with term of use.
- U.S. patent No.7. 198. 255
EU, Korea & Japan : Patent pending.
- Contents per box : 10 bags

Calm Pure Humidifier System (Oxygen therapy)



- Soft bag will save the storage space and reduce the garbage volume.



- Upper face of the adapter of the bag is made of soft film, it is easy to stick and connect with adapters.



- Uniform bubble can be acquired by precisely processed bubble generation part(bottom of the bag), efficient and quiet humidification can be provided.

Calm Pure Cold Nebulizer System (Aerosol therapy)



- The various particles size of 0.3~8micron can humidify evenly from the bronchial tree to the alveoli
- Oxygen level of inhalation can be set up at 35%, 40%, 45%, 50%, 60%, 80% and 98%.
- Large particle of water will be returned to the bag through recirculation tube.
- Minimum oxygen flow rate is 4L/min.
- Warm nebulization can be provided with optional aerosol heater.